

ARC® 2.0 SHOULDER BRACE

HCPCS: L3960

PRODUCT DESCRIPTION:

Our signature aluminum waistband is moldable to each patient's unique torso shape and prevents anterior migration that is common with shoulder braces. The ARC® 2.0 universal sling design folds to fit every patient with one brace - right or left, football player or gymnast. The ARC® 2.0 features 2.0 material for maximum breathability and comfort. 2.0 material moves moisture away from the skin to another layer of fabric for quick evaporation.



COMMON EXAMPLES OF USE:

- For non-surgical or post-op shoulder conditions that require diagnosis-specific arm positioning
- Rotator cuff repairs
- Bankart lesions
- SLAP lesions
- Glenohumeral dislocations / subluxation
- Posterior / Anterior capsule repairs
- Shoulder instabilities
- Muscle and tendon repair
- Joint reconstruction

PRODUCT FEATURES:

- Brace positions include gunslinger, neutral plane and external rotation
- Universal sling design to fit every patient with one brace
- Unique pistol grip adjusts with quick-pull tabs and keeps the arm from migrating forward out of the sling
- Optional underarm strap relieves pressure on the neck for larger patients and patients positioned in external rotation
- One-hand buckles, ease patient reapplication
- Cold therapy cutout in the sling

ORTHOMED, LLC
(323) 869-0000
info@orthomed-llc.com
9925 Painter Ave. Suite #D Whittier, CA 90605

ORTHO MED
Durable Medical Equipment