

# TRI-WAVE™ ELECTROTHERAPY STIMULATOR

HCPCS: E1399



## BENEFITS OF ELECTRIC STIMULATION:

- Pain Management Control
- Relax Muscle Spasms
- Increase Blood Circulation
- Reduce Inflammation
- Reduce Chronic or Post-Op Pain
- Re-Educate and Strengthen Muscles
- Increase Range of Motion (ROM)
- Prevent or Retard Disuse Atrophy

## PRODUCT DESCRIPTION:

The Tri-Wave™ Electric Stimulator is a multi-modality prescription only device with Interferential, TENS, and NMES treatment modalities used as an approach to pain management control, deep tissue penetration, and functional restoration. The Tri-Wave™ is designed to provide a portable user-friendly experience for the patient, while offering a myriad of possibilities for the clinician.

The CA ODG Guidelines supports patient use when there is evidence of increased functional improvement, less reported pain and evidence of medication reduction. ACOEM Guidelines, Chapter 3 states that “Electrical Stimulation can keep symptoms at bay temporarily, diminishing pain long enough so that patients begin to mobilize”.

## Electrotherapy and Heat Therapy Treatment

The Tri-Wave™ Electric Stimulator combined with the Thermotex Heating Pad is the perfect combination for pain management control. Combining electrotherapy and heat therapy modalities will keep the patients active and increase functional improvement.



## FEATURES:

- Two-channel protocols (4-lead output)
- Built-in rechargeable battery
- AC power adapter with built-in charger & indicator
- Five preset protocols including IF, TENS and NMS
- Optional programmable custom clinician protocols
- User-friendly menu interface
- Countdown treatment timer display
- Clinical compliance metering built-in
- Electrode leads have TouchProof\* output connectors

Tel: (888) 346-0888 • Fax: (866) 677-4106  
Email: info@orthomed-llc.com | www.orthomed-llc.com

**ORTHO MED**  
Durable Medical Equipment