

THIS DEVICE IS **TSA COMPLIANT**

This card should accompany your Ortho THERM FAR Infrared Heat Pad™
when traveling on an airplane

OrthoTHERM Smart Heating Pad™ is an electric
heating pad that uses **non-removable** lithium-ion
batteries that are fully compliant with TSA regulations.

Battery Specifications

Voltage: 14.4V

Charge: 5000mAh

Energy: 72Wh

NRTL Certified

SCAN ME

for Battery
Material Safety
Data Sheet
(MSDS)



This device can be taken onboard or checked on a commercial aircraft,
ensuring hassle-free travel.

Please visit **www.tsa.gov** for more information

Updated: May 2025

www.orthomed-llc.com